

BED-MON

Congratulations and thank you for choosing the BED-MON series of devices, your new BED-MON has been engineered to reduce false alarms due to movements in the bed to a minimum. These instructions have been written to help you with the basic operation of all models.

The BED-MON Plus model has additional features such as bed lamp function, and has a second output for those special requirements. The additional functions of BED-MON Plus will be covered separately in an included leaflet with BED-MON Plus.

The output plug has not been included in the pictures below as it is dependant on your brand and model of nurse call system installed in your complex, dependant on which model has been ordered a radio transmitter maybe installed inside of your new BED-MON.



BED-MON P (pendant jack)

Standard BED-MON

Smartlink Version with radio transmitter

INITIAL SETUP

ALL BED-MON models.

1. Insert the bed mat plug into the large BED-MON socket if you have BED-MON with pendant/floor mat socket then insert into the pendant socket.
2. Insert the power supply into 240 Volt power point and the other end into the small socket on BED-MON.
3. Either insert the corded plug that suits your nurse call system into the pendant socket of your nurse call wall plate **OR** follow the radio section below to setup a BED-MON with a radio option.
4. Turn the BED-MON on via the power point switch.

BED-MON with Radio Transmitter.

First time use you will need to program the BED-MON radio transmitter into a receiving device such as medical dialler etc. Follow the instruction as per your medical dialler or similar on the procedure regarding how to program the BED-MON radio transmitter into your medical dialler or similar.

Use the programming button on the top of the BED-MON to add the BED-MON into your radio receiving device, the same button can be used to test the transmitter once programming has been completed.

Goto step 4.



Blanking Plug

OPERATION

Your new BED-MON does not need resetting as it resets itself when ever a person has returned to bed. However, you will still need to follow your normal operating procedures regarding the cancelling of any Nurse call system alarms.

The BED-MON only generates a pulse to activate your nurse call system and will not clear or cancel any call once the BED-MON has activated your nurse call system.

It is recommended that the BED-MON not be turned on until the person is in bed. Once the person has settled into bed then it is OK to turn on the BED-MON.

PLEASE NOTE:

1. Your BED-MON may generate an alarm when first switched on and has not been fully setup with a person in bed and all cables plugged in.
2. With no bed mat plugged in the BED-MON will generate an alarm (this is normal) unless an optional blanking plug is used instead of the bed mat.
3. If you do not want to turn off the BED-MON it is mandatory that you use a blanking plug to stop any alarms.
4. The pendant cord will still function if you do turn off the BED-MON. Also the pendant socket does not need a blanking plug when pendant cord not in use.

U4 2800 Logan Road, Underwood, Qld 4127
PO Box 1262, Springwood, Qld 4127
www.bfelectronics.com.au
sales@bfelectronics.com.au
Office: +61 7 3117 0552
Mobile: +61 478 068 088



BFELECTRONICS